English**Description**

Prostrolane Blanc_B is a transparent gel supplied in a glass syringe. The product is for single use only. Prostrolane Blanc_B is a sterile medical device, dermal resorbable implant that contains a non-cross-linked Sodium Hyaluronate and peptide complex (Oligopeptide-92, Octapeptide-11, Feruloyl Oligopeptide-33, Valprooyl Oligopeptide-33) Sodium Hyaluronate is non-animal origin and biocompatibility material. The role of Sodium Hyaluronate in the skin is to deliver peptide, to hydrate the skin by holding in water and to act as a cushioning agent. Peptide complex consists of 4 kinds peptides that short chain of amino acids linked by peptide bonds. The role of the peptide complex is to enhance the moisturizing effect, to improve the appearance of photo-aged skin and to delay the degradation of Sodium Hyaluronate. The sodium hyaluronate peptide complex is degradable within 4 weeks. Each carton contains 1 prefilled syringe of 2ml Prostrolane Blanc_B and two traceability labels. (one to be given to the patient, and another to be kept by the doctor in the patient's file)

Composition

PBS/Phosphate Buffered Saline), Sodium Hyaluronate (0.7%), Peptide Complex (Oligopeptide-92, Octapeptide-11, Feruloyl Oligo peptide-33, Valprooyl Oligopeptide-33)

Intended purpose

Prostrolane Blanc_B is indicated for dermal implantation or the treatment of damaged subcutaneous tissue and wrinkles caused by chemical, physical and traumatic effects, as well as sun damaged tissue.

- Site of Application: Face, neck

- Recommended usage volume : 2ml / area (Maximum: 4ml / area)

- Frequency : 1 session / 2 weeks (1 protocol = 2-4 sessions average) (repeat protocol as necessary)

Treatment procedure

Pretreatment Guidelines
Prior to treatment, the patient should avoid taking aspirin, nonsteroidal anti-inflammatory medications, St. John's Wort, or high doses of Vitamin E supplements. These agents may increase bruising and bleeding at the treatment site.

1. It is necessary to counsel the patient and discuss the appropriate indication, risks, benefits and expected responses to the Prostrolane Blanc_B treatment. Advise the patient of the necessary precautions before commencing the procedure.

2. Assess the patient's need for appropriate anesthetic treatment for managing comfort, i.e., topical anesthetic, local or nerve block.

3. The patient's face should be washed with soap and water and dried with a clean towel.

Cleanse the area to be treated with alcohol or another suitable antiseptic solution.

4. Sterile gloves are recommended while using Prostrolane Blanc_B.

5. Before treatment, press rod carefully until a small droplet is visible at the tip.

6. Prostrolane Blanc_B is administered using a thin gauge needle (27G x 13mm or 30G x 13mm) that is sterilized. The sterilized needle is inserted at an approximate angle of 30° parallel to the length of the injected site. Prostrolane Blanc_B should be inserted to the mid to deep dermis. If Prostrolane Blanc_B is injected too superficially this may result in visible lumps and/or bluish discoloration.

7. Inject Prostrolane Blanc_B applying even pressure on the plunger rod. It is important that the injection is stopped just before the needle is pulled out of the skin to prevent material from leaking out or ending up too superficially in the skin.

8. Only correct to 100% of the desired volume effect. Do not overcorrect. If cutaneous deformities the best results are obtained if the defect can be manually stretched to the point where it is eliminated. The degree and duration of the correction depend on the character of the defect treated.

9. Typical usage for each treatment session is specific to the site as well as wrinkle severity. Or to apply on the cleared skin surface, to wait for full absorption.

10. After treatment, stick the traceability label on the patient's medical record.

Direction for Assembly
Assembly by 27G or 30G needle(hypodermic needle) to syringe.

For safe use of Prostrolane Blanc_B, it is important that the needle is properly assembled. Improper assembly may result in separation of the needle and syringe during injection.


Storage condition
Store between 2°C - 25°C. Protect from freezing and sunlight. Make sure there are no visible signs of damage to the packaging before use.

Shef life
2 years under recommended conditions.

A Manufactured by Caregen Co., Ltd.
46-38 Caregen Bldg., LS-ro 91beon-gil,
Dongan-gu, Anyang-si, Gyeonggi-do 14119 South Korea
TEL +82-31-420-9200
Fax +82-31-452-3869
e-mail : sales@caregen.com

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Contraindications

Prostrolane Blanc_B must not be used:
- for the combination with peeling, laser treatment or ultrasound;
- if the patient has cutaneous disorder, inflammation or infection at the treatment site or near to this site;
- in the case of patients have a known hypersensitivity to hyaluronic acid, with a history of severe allergy or anaphylactic shock;
- in case of the patients with autoimmune diseases;
- Due to possible interactions with other filling implants, which have not been researched, it is inadvisable to use Prostrolane Blanc_B into sites in the presence of other filling implants;
- in pregnancy, breast-feeding mother, and in children or minors under 18 ages;
- for patients with bleeding disorders, skin color disorders

Composition

PBS(Phosphate Buffered Saline), Sodium Hyaluronate (0.7%), Peptide Complex (Oligopeptide-92, Octapeptide-11, Feruloyl Oligo peptide-33, Valprooyl Oligopeptide-33)

Intended purpose

Prostrolane Blanc_B is indicated for dermal implantation or the treatment of damaged subcutaneous tissue and wrinkles caused by chemical, physical and traumatic effects, as well as sun damaged tissue.

- Site of Application: Face, neck

- Recommended usage volume : 2ml / area (Maximum: 4ml / area)

- Frequency : 1 session / 2 weeks (1 protocol = 2-4 sessions average) (repeat protocol as necessary)

Treatment procedure

Pretreatment Guidelines
Prior to treatment, the patient should avoid taking aspirin, nonsteroidal anti-inflammatory medications, St. John's Wort, or high doses of Vitamin E supplements. These agents may increase bruising and bleeding at the treatment site.

1. It is necessary to counsel the patient and discuss the appropriate indication, risks, benefits and expected responses to the Prostrolane Blanc_B treatment. Advise the patient of the necessary precautions before commencing the procedure.

2. Assess the patient's need for appropriate anesthetic treatment for managing comfort, i.e., topical anesthetic, local or nerve block.

3. The patient's face should be washed with soap and water and dried with a clean towel.

Cleanse the area to be treated with alcohol or another suitable antiseptic solution.

4. Sterile gloves are recommended while using Prostrolane Blanc_B.

5. Before treatment, press rod carefully until a small droplet is visible at the tip.

6. Prostrolane Blanc_B is administered using a thin gauge needle (27G x 13mm or 30G x 13mm) that is sterilized. The sterilized needle is inserted at an approximate angle of 30° parallel to the length of the injected site. Prostrolane Blanc_B should be inserted to the mid to deep dermis. If Prostrolane Blanc_B is injected too superficially this may result in visible lumps and/or bluish discoloration.

7. Inject Prostrolane Blanc_B applying even pressure on the plunger rod. It is important that the injection is stopped just before the needle is pulled out of the skin to prevent material from leaking out or ending up too superficially in the skin.

8. Only correct to 100% of the desired volume effect. Do not overcorrect. If cutaneous deformities the best results are obtained if the defect can be manually stretched to the point where it is eliminated. The degree and duration of the correction depend on the character of the defect treated.

9. Typical usage for each treatment session is specific to the site as well as wrinkle severity. Or to apply on the cleared skin surface, to wait for full absorption.


10. After treatment, stick the traceability label on the patient's medical record.

Direction for Assembly
Assembly by 27G or 30G needle(hypodermic needle) to syringe.

For safe use of Prostrolane Blanc_B, it is important that the needle is properly assembled. Improper assembly may result in separation of the needle and syringe during injection.

Storage condition
Store between 2°C - 25°C. Protect from freezing and sunlight. Make sure there are no visible signs of damage to the packaging before use.

Shef life
2 years under recommended conditions.

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