DERMAFIX GENOSYS

GENE

DTS MG Co., Ltd.

GENOSYS DERMAFIX

"GENOSYS DERMAFIX is an automatic device for microneedling therapy."

It helps increase the absorption rate of active ingredients of skin care products by creating physical pathways through skin. And it helps collagen and elastin induction through natural wound healing process of skin.



GENOSYS DERMAFIX needle cartridge is designed for GENOSYS DERMAFIX auto microneedling therapy. It is sterilized for safety and special design of needle cartridge prevents blood or solution from contaminating the machine.



Features of GENOSYS DERMAFIX

•Disc type needle: Each metal disc has four needles. And four discs compose one needle cartridge. When pricking the skin with the disk needles, the needles do not give severe damage to the skin and can be prevented from being broken or fell out.

• Number of needle-16: The number of needle is 16, thus making more micro channels when compared to other autoneedling devices whose needle cartridges have 9-12 needles in general.

- Thickness of needle- 0.20mm: The pricking sensation is minimized.
- Cutting shape of needle-hexagon:

The circumferential surface of the drum is configured to form a polygonal cross-section of at least six sides instead of a circular cross section, so that, when pricking the skin with the needles, the contact surface of each needle relative to the skin can be reduced, thus minimizing the pain.

Hexagon

- Sterilized needles: hygiene, safety
- Adjustable motor power and speed: 9 levels

• Adjustable needle length-0.25mm~2.00mm: According to the purpose of treatment, needle length can be adjusted, which allows for more effective and affordable treatment for a patient with a multitude of skin conditions requiring different length needles.→ cost-effective

Effects of GENOSYS DERMAFIX

- Skin firming and lifting
- Skin brightening and whitening
- Improvement of fine lines and wrinkle
- Improvement of acne scar
- Improvement of burn scar
- Improvement of stretch mark
- Treatment of hair loss (alopecia)
- Contraction of pores

Specifications of GENOSYS DERMAFIX

- Composition: device, adapter, needle cartridge (2 pcs)
- RPM: 1000~6000
- Adapter: Input 100V~240V (50/60Hz) Output DC 1~12V (100~1000mA)
- Motor voltage: 12V Weight of device: 101g









Capactuel L'esthétique Préventive & Réparatrice

Tél.: 01 85 34 31 37 contact@capactuel.com

www.capactuel.com